


veterinary technology for life

# Eickemeyer®

[www.eickemeyer.hu](http://www.eickemeyer.hu)

Telephone 06 1 206 9555

## *Experience Endoscopy*



# EICKVIEW SD ENDOSCOPY SYSTEM

## Unlock the future of veterinary care with the EickView SD Endoscopy System!

This state-of-the-art system features a high-quality endoscopy monitor with a large, intuitive touchscreen, seamlessly connecting to three versatile “semi-disposable” endoscopes. This system allows the veterinary practitioner to experience the benefits of flexible endoscopy (Gastroscopy, Bronchoscopy, Rhinoscopy, etc.) without the hefty investment of traditional video endoscopes. Each endoscope is equipped with a high-resolution CMOS image sensor and an LED light source, ensuring crystal-clear visuals.

### EickView SD Endoscopy Monitor

- ▶ 15.6" touchscreen with 1,920 x 1,080 Full HD resolution
- ▶ Functions: White balance, zoom, freeze, photo (JPG), video recording (MP4) and DICOM
- ▶ 64 GB storage capacity
- ▶ Mains and battery operation, battery life up to 5 hours
- ▶ Outputs: HDMI, DVI, SDI and USB
- ▶ Dimensions (in cm): W 39 x H 25 x D 4
- ▶ Weight: 2.7 kg

**307000**

### EickView 150SD

- ▶ Working length: 150 cm
- ▶ Ø outside: 6.0 mm
- ▶ Ø working channel: 2.8 mm
- ▶ Direction of view: 0°
- ▶ Viewing angle: 120°
- ▶ Bending up / down: 180°
- ▶ Bending right / left: 120°

**307010**

### EickView 60SD

- ▶ Working length: 60 cm
- ▶ Ø outside: 3.2 mm
- ▶ Ø working channel: 1.2 mm
- ▶ Direction of view: 0°
- ▶ Viewing angle: 120°
- ▶ Bending up / down: 210°

**307020**

### EickView 60SD-XS

- ▶ Working length: 60 cm
- ▶ Ø outside: 2.7 mm
- ▶ Ø working channel: 1.2 mm
- ▶ Direction of view: 0°
- ▶ Viewing angle: 120°
- ▶ Bending up / down: 210°

**307030**

### EickView 65SD-URS

- ▶ Working length: 65 cm
- ▶ Ø outside: 2.1 mm
- ▶ Ø working channel: 1.2 mm
- ▶ Direction of view: 0°
- ▶ Viewing angle: 120°
- ▶ Bending up / down: 270°

**307040**



## EICKVIEW 5-IN-1 ENDOSCOPY SYSTEM

**ORDER  
ONLINE**

www.eickemeyer.hu

### EickView 5-in-1 Endoscopy System

The EickView 5-in-1 endoscopy system is a true all-rounder for veterinary medicine. The easy-to-use, mobile device offers the option to work with both, flexible and rigid endoscopes. This means that all types of endoscopy can be covered with just one device – from diagnosis to surgical procedures.

The EickView 5-in-1 combines a large, 15.6" Full HD display, a powerful LED light source, a connection for flexible endoscopes from the EickView HD series with air and water pump as well as digital data storage via SD Card in one system.

**306200**

### Outdoor case for EickView 5-in-1

- ▶ With foam insert
- ▶ Made of robust plastic

**30620001**

### EickView HD and HD+ Video Endoscopes

The EickView HD and HD+ Video Endoscopes deliver brilliant images for bronchoscopy, gastroscopy and various other applications. The new HD+ Endoscopes are characterised by their reduced outer diameters and enlarged working channels, making them even more variable. Due to the large number of available scope sizes, they are suitable for all applications in small animals and equine.

### Technical specifications

- ▶ Viewing direction: 0°
- ▶ Distal end contains 1 CMOS Camera Chip, 2 Light Guides, Working Channel, Lens Rinsing Nozzle
- ▶ Scope of delivery: EickView Video Endoscope, Pressure Tester, 2 Cleaning Brushes, Cleaning Set, Carrying Case

### EickView 5-in-1 Camera Kit

The EickView 5-in-1 camera with zoom adapter and light guide cable completes the comprehensive package. It allows all standard rigid endoscopes to be connected to the main system.

**306201**



306201



306200

306268

### EickView HD Video Endoscopes

Item No.	Description	Working length	Outer Ø	Working channel	Depth of field	Field of view	Deflection
<b>306177</b>	EickView 60 HD Endoscope	60 cm	3.5 mm	1.2 mm	3 – 50 mm	120°	up: 160° / down: 130° right / left: 0°
<b>306178</b>	EickView 100 HD Endoscope	100 cm	6.0 mm	2.0 mm	3 – 50 mm	120°	up: 210° / down: 90° right / left: 100°
<b>306266</b>	EickView 150-6 HD+ Endoscope	150 cm	6.0 mm	2.8 mm	2 – 100 mm	120°	up: 210° / down: 90° right / left: 100°
<b>306268</b>	EickView 150-8 HD+ Endoscope	150 cm	8.0 mm	2.8 mm	3 – 50 mm	120°	up: 210° / down: 90° right / left: 100°
<b>306270</b>	EickView 350 HD+ Endoscope	350 cm	13.0 mm	3.7 mm	3 – 100 mm	140°	up: 180° / down: 180° right / left: 160°

## EICKVIEW PORTABLE LED VIDEO ENDOSCOPES

Due to their extremely compact dimensions, the EickView LED Video Endoscopes are the perfect solution for broncho-, gastro-, rhino-, tracheo-, uretero- and vaginoscopy in dogs and cats and for the examination of the upper respiratory tract in horses. The devices work independently, as the powerful LED light source and the high resolution camera of the stand-alone systems are integrated in the distal end of the endoscopes. In addition, the folding 3.5" monitor allows both, direct video and image playback as well as storage via SD card. Furthermore the EickView LED Video Endoscopes can be connected to a separate monitor.

### Technical specifications

- ▶ Folding 3.5" monitor for video and picture playback
- ▶ Recording and archiving on integrated SD card
- ▶ File formats: JPG, AVI
- ▶ Integrated powerful LED light source
- ▶ Rechargeable Battery with up to 3 hours of operation
- ▶ USB port and video out
- ▶ Depth of field: 5 – 80 mm
- ▶ Direction of view: 0°
- ▶ Field of view: 110°



306166

EickView portable LED Video Endoscopes					
Item No.	Description	Working length	Outer Ø	Working channel	Deflection
<b>306009</b>	EickView 60S Endoskope	60 cm	3.3 mm	1.2 mm	up 180° down 130°
<b>306011</b>	EickView 60 Endoskope	60 cm	4.2 mm	1.3 mm	up 180° down 130°
<b>306010</b>	EickView 70 Endoskope	70 cm	6.0 mm	2.0 mm	up 180° down 130°
<b>306166</b>	EickView 100E Endoskope	100 cm	8.0 mm	2.8 mm	up 150° down 130°
<b>306164</b>	EickView 150E Endoskope	150 cm	8.0 mm	2.8 mm	up 130° down 130°

## ACCESSORIES EICKVIEW ENDOSCOPES


		EickView 60 / 60S/60 HD / 60SD / 60SD-XS	EickView 70	EickView 100 HD/150 HD	EickView 100E/150E/150L / 150-8 HD+ / 150-6 HD+ / 150SD	EickView 350 HD+
<b>Biopsy Forceps</b>	Item No.	<b>306311</b>	<b>306314</b>	<b>306318</b>	<b>306316</b>	<b>306315</b>
	Length (cm)	120	120	230	230	400
	Ø (mm)	1.0	1.8	1.8	2.3	2.2
<b>Alligator Forceps</b>	Item No.	<b>306211</b>	<b>306214</b>	<b>306217</b>	<b>306216</b>	<b>306215</b>
	Length (cm)	120	120	230	230	400
	Ø (mm)	1.0	1.8	1.8	2.4	2.2
<b>Grasping Forceps</b>	Item No.	<b>306340</b>	<b>306346</b>	<b>306342</b>	<b>306343</b>	<b>306345</b>
	Length (cm)	120	120	230	230	400
	Ø (mm)	1.0	1.8	1.8	2.3	2.2
<b>Rat Tooth Forceps</b>	Item No.	<b>306431</b>	<b>306435</b>	<b>306434</b>	<b>306432</b>	<b>306437</b>
	Length (cm)	120	120	230	230	400
	Ø (mm)	1.0	1.8	1.8	2.3	2.2
<b>Basket Retrieval Device</b>	Item No.	<b>306336</b>	<b>306332</b>	<b>306331</b>	<b>306329</b>	<b>306335</b>
	Length (cm)	120	120	230	230	400
	Ø (mm)	1.0	2.0	2.0	2.3	2.45
	Basket size (mm)	20	20	40	40	50
	No. of filaments	4	6	6	6	6
<b>Polypectomy Snare</b>	Item No.		<b>306426</b>	<b>306422</b>	<b>306423</b>	<b>306421</b>
	Length (cm)		120	230	230	400
	Ø (mm)		1.8	1.8	2.4	2.45
<b>Cytology Brush</b>	Item No.	<b>306415</b>	<b>306416</b>	<b>306417</b>	<b>306418</b>	<b>306414</b>
	Length (cm)	120	160	230	230	400
	Ø (mm)	1.0	1.8	1.8	2.6	2.5
<b>Cleaning Brush</b>	Item No.	<b>306444</b>	<b>306440</b>	<b>306441</b>	<b>306442</b>	<b>306445</b>
	Length (cm)	120	160	230	230	400
	Ø (mm)	2.2	3.0	3.0	5.0	7.0
<b>Suction Catheter</b>	Item No.	<b>306353</b>	<b>306353</b>	<b>306352</b>	<b>306352</b>	<b>306355</b>
	Length (cm)	100	100	350	350	400
	Ø (mm)	1.1	1.1	1.8	1.8	2.45



306311




306211




306340



306431




306336




306426



306415



306444



306353



# LAPAROSCOPY BASIC SET


The Laparoscopy Basic Set provides you with a solid foundation for minimally invasive procedures. It enables gentle treatments in areas such as the adrenal glands, gallbladder, lungs, intestines, mediastinum, and reproductive tract.

The set can be easily expanded with specialized instruments and accessories to meet the individual needs of your practice. This keeps you flexible and optimally positioned to perform your surgical procedures efficiently and precisely.


## Complete set, consisting of:

- ▶ Rigid Endoscope, L 32.8 cm, Ø 5 mm, 30° (Item No. 305004)
- ▶ EickView 5-in-1 Endoscopy System (Item No. 306200)
- ▶ EickView Camera Kit (Item No. 306201)
- ▶ Trocar Sleeve, L 9.5 cm, inner Ø 5.5 mm (Item No. 30495000)
- ▶ Safety Trocar, L 16 cm, Ø 5.5 mm (Item No. 30495031)
- ▶ Trocar Sleeve, L 10.1 cm, inner Ø 11 mm (Item No. 30500005)
- ▶ Safety Trocar, L 17.1 cm, Ø 11 mm (Item No. 30495033)
- ▶ Veress cannula, L 10 cm, Ø 2.0 mm (Item No. 305002)
- ▶ EickTron 160 W electrosurgical unit
  - Output power adjustable via foot or hand switch
  - Max. Output: 160 W (cut), 100 W (coagulate), 60 W (bipolar)
 (Item No. 323135)
- ▶ Bipolar connecting cable (Item No. 323106)
- ▶ Bipolar connecting adapter (Item No. 323139)
- ▶ EickCut bipolar coagulation grasping and cutting forceps, L 32 cm, Ø 5 mm (Item No. 30495013)
- ▶ Blades for EickCut forceps, 5 pieces per pack (Item No. 30495014)
- ▶ HF connecting cable for monopolar instruments (Item No. 30495050)
- ▶ Monopolar universal grasping forceps L 36 cm, Ø 5.0 mm (Item No. 30495016)
- ▶ CO<sub>2</sub> insufflator
  - Pressure Range: 1 – 30 mmHg
  - Five Insufflation Modes
 (Item No. 305300)
- ▶ Pressure reducer for CO<sub>2</sub> insufflator (Item No. 30526001)
- ▶ Cannula Snap Reducer Adapter, Ø 10 or 11 mm to 5,5 mm (Item No. 30500007)


## 306040



305004



30495013



30495033



30500005



30495000



323106


30495031

30495033

305002



323139



30495050



30495014



30495016



305300



30526001



306200



306201



306201



306201

## EICKEMEYER® LAPAROSCOPY POSITIONER

**ORDER  
ONLINE**  
www.eickemeyer.hu

The EICKEMEYER® Laparoscopy Positioner enables greater visualisation as the patient can be rotated. The rotation helps major abdominal organs to fall away from the body wall, leaving the ovary in plain view and ready for removal.


### Features

- ▶ Tilt and secure patient in different steps both left and right
- ▶ Easy to assemble and disassemble
- ▶ Easy to clean

### EICKEMEYER® Laparoscopy Positioner (Dogs)

- ▶ Dimensions (in cm): L 137 x W 44 x H 34,5

**305100**



305100

### EICKEMEYER® Laparoscopy Positioner (Cats)

- ▶ Features ability to tilt at a greater angle for small patients
- ▶ Dimensions (in cm): L 73 x W 28.5 x H 23

**305110**



305110

## ENDOSCOPY TROLLEYS

### Endoscopy Trolley Stainless steel

- ▶ Three adjustable shelves: 54 x 50 cm
- ▶ One fixed bottom shelf: 41 x 57 cm
- ▶ Border: 1.5 cm
- ▶ Four lockable castors
- ▶ Dimensions (in cm): H 170 x W 60 x D 60

**306160**



306160

### Endoscopy Trolley

- ▶ Four storage shelves 60 x 49 cm + one drawer
- ▶ Raised edge: 10 mm
- ▶ Two energy columns for power supply
- ▶ Standard color: white / grey
- ▶ Four conductive twin castors
- ▶ Dimensions (in cm): H 139.5 x W 60 x D 62

**306150**



306150



306155

### EICKEMEYER® Endoscopy Trolley

- ▶ Three adjustable shelves: 54 x 50 cm
- ▶ One fixed bottom shelf: 50 x 54 cm
- ▶ LCD monitor shelf: 30 x 50 cm
- ▶ Standard colour: green / grey
- ▶ Four lockable castors
- ▶ Dimensions (in cm): H 170 x B 60 x T 70

**306155**

## EICKVIEW COMPACT ENDOSCOPE SET

The specially developed EickView Compact Endoscope provides the perfect solution for rhinoscopy, tracheoscopy and vaginoscopy in small animals due to its small diameter (3.8 mm) and integrated working channel.

- ▶ Highly mobile
- ▶ Smooth and minimally invasive distal end
- ▶ Integrated working channel Ø 1.3 mm
- ▶ Outer diameter: 3.8 mm
- ▶ Working length 14 cm
- ▶ Viewing angle: 30°
- ▶ Opening angle 90° (wide-angle)
- ▶ Foreign body forceps Ø 1 mm, L 25 cm
- ▶ Attachable LED Light Source
- ▶ Long battery life and spare battery with an operating time of 4 hours each
- ▶ Three Adjustable light levels
- ▶ Stopcocks: two integrated stopcocks allow steady insufflation and irrigation
- ▶ Connection: STORZ, WOLF or ACMI connectors
- ▶ Eyepiece: DIN ocular allows the connection of all standard C-mount cameras
- ▶ Set includes: EickView Compact endoscope incl. fitted otoscopy protection aid, portable LED light source incl. charger, foreign body forceps, spare battery, 2 cleaning brushes, padded carrying case

**306122**



## EICKVIEW ULTRA COMPACT ENDOSCOPE SET

The EickView Ultra Compact is the new generation of multifunctional rigid endoscopes. It is ideal for all ENT applications as well as tracheo- and vaginoscopies with its 3.0 mm diameter.

- ▶ Smooth and minimally invasive distal end by fusion of endoscope and sheath
- ▶ No dirt deposits in the distal end
- ▶ Small outer diameter allows easy access
- ▶ Integrated working channel enables the removal of foreign bodies and biopsies
- ▶ Integrated working channel Ø 1.3 mm
- ▶ Outer diameter: 3.0 mm
- ▶ Working length: 11 cm
- ▶ Viewing angle: 0°
- ▶ Opening angle: 90° (wide-angle)
- ▶ Endoscope can be cold sterilised
- ▶ Stopcocks: two integrated stopcocks allow steady insufflation and irrigation
- ▶ Set includes: EickView Ultra Compact endoscope, portable LED light source incl. charger, foreign body forceps, spare battery, 2 cleaning brushes, padded carrying case

**306112**

